

**NEW:**

# derma Litecheck® UV Multimedia

now also available with integrated USB camera

**EUR 459,00**

plus shipping costs outside Germany and legal VAT

- **Operation: e.g. with Windows 10 (Camera and Media Player)**
- **Editing of films (cut etc.) with common Windows-programmes**

Multimedia device with integrated USB camera,  
encl. connecting cable CMOS, 3.14 Megapixel, 80° V5 Fixfokus, USB 2.0

## **Initiating operation:**

- Insert mains plug into the socket, which is located at the left side of the device. Switch on (switch lights up red).
- Remove protective cap from camera eye, which is integrated in the glass plate under the electronic.
- Insert the USB cable into the corresponding socket on the back side of the device.
- Connect the USB cable to your computer
- Sound Signal/screen note: „USB camera recognized/installed“
- Example Windows 10: open programme „Camera“

**NEW:**

# derma Litecheck® UV Multimedia

now also available with integrated USB camera

## visualize and communicate

by means of the smart desk solution made by Dermalux®

### **Examples for application:**

scientific studies  
design of lectures  
animation with pictures and/or videos  
live-switch via Skype  
presentations by laptop/beamer/monitor/TV  
large projections

personnel training  
trade fair  
congress  
group training  
prove of evidence  
customer presentation  
sales meetings

archiving of data  
long-term investigations

and much more...

